



Myth
VS. **FACT**

Exposing the Myths of ACL Reconstruction


ARTHROTEK®
A BIOMET COMPANY

Biomechanical Myths


Myth #1

Bone-patellar tendon-bone is the strongest and stiffest graft.

Fact: Bone-patellar tendon-bone (BPTB) is one-half the strength of double-looped semitendinosus and gracilis (DLSTG).^{1,3,6,8}



Fact: BPTB is half as stiff as DLSTG.^{1,3,6,8}



Conclusion: The literature concludes that the double-looped semitendinosus and gracilis is the strongest and stiffest graft.^{1,3,6,8}

Myth #2

Interference screw fixation of a bone patellar tendon-bone (BPTB) graft is the "gold standard" for ACL reconstruction.

Fact: The strength of BPTB fixed with interference screws is 416N.¹

Fact: The stiffness of BPTB fixed with interference screws is 51N/mm.¹

Fact: Slippage of BPTB graft fixed with interference screws is an excessive 3.8mm at 500N.²

Fact: DLSTG fixed with Bone Mulch™ Screw and WasherLoc™ strength is 903N.⁷

Fact: DLSTG fixed with Bone Mulch™ Screw and WasherLoc™ stiffness is 208N/mm.

Fact: DLSTG fixed with Bone Mulch™ Screw and WasherLoc™ slippage is .8mm in porcine bone. DLSTG fixed with Bone Mulch™ Screw and WasherLoc™ slippage is 1.8mm in human bone.³

Bone Mulch Screw vs. Interference Screw for Femoral Fixation^{5,7}

| Type of Fixation | Yield Load | Stiffness |
|---------------------------|------------|-----------|
| BPTB w/Interference Screw | 416N | 51N/mm |
| DLSTG w/Bone Mulch Screw | 1125N | 225N/mm |

WasherLoc vs. Interference Screw for Tibial Fixation³

| Type of Fixation | Yield Load | Stiffness |
|---------------------------|------------|-----------|
| BPTB w/Interference Screw | 340N | 248N/mm |
| DLSTG w/WasherLoc | 905N | 273N/mm |

Conclusion: The literature concludes that interference screw fixation is not the gold standard for ACL reconstructions.

Rehabilitation Myths

Myth #3

Aggressive rehabilitation can be applied with all fixation methods.

Fact: Graft load can be as high as 560N with a 1500N quadriceps contraction.²

Fact: Endobutton femoral suspended from a tape fails at 430N.⁸

Fact: Endobutton femoral stiffness is 23N/mm.⁸

Fact: Interference screws in the tibia fails at 340N.⁴

Fact: The stiffness of an interference screw in the tibia is 248N/mm.⁶

Fact: Interference screws in the tibia slip 3.7mm at 500N.⁶

Fact: Tibial fixation with a 20mm Washer fails at 724N.⁶

Fact: Tibial fixation with a 20mm Washer stiffness is 126N/mm.⁶

Fact: Tibial fixation with a 20mm Washer slips 3.5mm at 500N.⁶

Fact: Belt-buckle/staples slips 3.3mm at 500N.⁶

Fact: Belt-buckle/staples on the tibia stiffness is 118N/mm.⁶

Fact: Belt-buckle/staples slips 3.3mm at 500N.⁶

Conclusion: The literature concludes that aggressive rehabilitation of BTB cannot be applied with all fixation methods.

Myth #4

Aggressive rehabilitation should not be used with hamstring grafts.

Fact: The DLSTG is two times stronger than BPTB grafts.^{1,3,6,8}

Fact: DLSTG is stiffer than BPTB grafts.^{1,3,6,8}

Fact: The Bone Mulch™ Screw/WasherLoc™ graft fixation complex meets all criteria to aggressively rehabilitate.

Fact: Extension exercises can induce loads in ACL greater than strength of most fixation methods.

Conclusion: The literature concludes that aggressive rehabilitation of hamstrings can be effective following the use of Bone Mulch™ Screw and WasherLoc™ fixation of the DLSTG.

References

- ¹ Brown, C.H., Jr.; Steiner, M.E., and Carson, E.W.: The use of hamstring tendons for anterior cruciate ligament reconstruction. Technique and results. *Clinical Sports Medicine*, 12: 723-56, 1993.
- ² Goss, Howell, Hull: Quadriceps Load Aggravates and Roofplasty Mitigates Active Impingement of Anterior Cruciate Ligament Grafts Against the Intercondylar Roof. *Journal of Orthopaedic Research*, 1998.
- ³ Howell, S.M., and Hull, M.L.: Aggressive rehabilitation using hamstring tendons: Graft construct, tibial tunnel placement, fixation properties, and clinical outcome. *American Journal of Knee Surgery*, 1998.
- ⁴ Howell, S.M. and Hull, M.L.: Structural properties of six tibial fixation methods for anterior cruciate ligament soft tissue grafts. *American Journal of Sports Medicine*, 27: 35-43, 1999.
- ⁵ Liu, S. H.; Kabo, J.M., and Osti, L.: Biomechanics of two types of bone- tendon- bone graft for ACL reconstruction. *Journal of Bone and Joint Surgery*, 77-B: 232-5, 1995.
- ⁶ Magen, H. E.; Howell, S.M., and Hull, M.L., Slippage and stiffness of six fixation devices used to secure a double- looped tendon soft tissue graft to the tibia, AAOS Annual Meeting, New Orleans, LA, 1998.
- ⁷ Rowden, N.J.; Sher, D.; Rogers, G.J., and Schindhelm, G.: Anterior cruciate ligament graft fixation: Initial comparison of patellar tendon and semitendinosus autografts in young fresh cadavers. *American Journal of Sports Medicine*, 25: 472- 478, 1997.
- ⁸ To, J. T.; Howell, S.M., and Hull, M.L.: Factors affecting the stiffness of anterior cruciate ligament replacements at implantation, AAOS Annual Meeting, Atlanta, GA, 1997.

Bone Mulch® Screw and WasherLoc™ were developed in conjunction with Stephen M. Howell, M.D., Sacramento, California.

WasherLoc™ and Bone Mulch™ Screw are trademarks of Arthrotek, Inc.

Arthrotek® is a registered trademark of Arthrotek, Inc.

ARTHROTEK®
A BIOMET COMPANY

P.O. Box 587, Warsaw, IN 46581-0587 • 219.267.6639 • ©2000 Arthrotek, Inc. All Rights Reserved
web site: <http://www.arthrotek.com> • eMail: arthrotek@arthrotek.com
Form No. Y-BMT-625R/101500/Intranet Only